

Staining Cells in a 96-well Plate (PBMCs)

1. Prepare a single cell suspension of up to 1×10^7 cells per mL.
2. Process cells with FC block prior to aliquoting out
 - a. For Human use FC Block (available from many vendors)
 - b. For mouse use CD16/CD32 purified.
3. Add 50uL of cells to each well.
4. Wash the cells by adding 200uL/well of wash buffer to each well.
5. Centrifuge 400-600g for 5 minutes.
6. Add a viability dye and incubate according to that dye's instructions
7. Wash the cells by adding 200uL/well of wash buffer to each well.
8. Centrifuge 400-600g for 5 minutes.
9. Add 50uL of correctly diluted Ab to each of your control wells.
10. Add 50uL of master-mix to each sample well.
11. Mix gently for 2-3 minutes on a plate shaker.
12. Incubate for 20-30 minutes on ice or at 2-8 degrees in the dark.
13. Wash the cells by adding 200uL/well of wash buffer to each well.
14. Centrifuge 400-600g for 5 minutes.
15. Discard supernatant.
16. Repeat wash – Steps 13-15
17. Resuspend cells in 100uL staining buffer and 100uL of fixation buffer.
18. Analyze samples as soon as possible but no later than 3 days.

Fixation Buffers:

BD – Cytotfix 554655

Biolegend – Fixation Buffer 420801

Thermofisher – Fixation buffer 08-8222-49

FC Block:

BD – Human FC Block 564219

BD – Mouse FC Block 553141S

Biolegend – Human FC Block 422302

Biolegend – Mouse FC Block 101320

Thermofisher – Human FC Block 14-9161-73

Thermofisher – Mouse FC Block 14-0161-82