

Core staff pre-training checklist:

- Assign project to training session in PPMS
 - Change user rights to autonomous for both instruments
 - Add lab login to Diva software (on both instruments) if none exists
 - o Set to automatically export FCS files after recording in User Preferences -> Export
 - Request HSRB badge access if needed
-

Two instruments (same configuration) available: Winship C-5027 and HSRB E-362

- Separate badge access request if using the A3 in Winship (instructions on ECFCC website)
 - o <https://sites.winship.emory.edu/facilities>
- First user of the day/instrument off: power on UV, warm-up period of 20 minutes is required
- Only BD falcon tubes may be used (VWR #352054)

Startup/Fluidics

- PPMS login
- Windows (BDAdmin) login
- Startup if first user/instrument off:
 - o Green button on right side of instrument
 - o UV laser (355nm) must be turned on and set to 60mW
 - 20min warmup period
 - Minimize Coherent connection
- Check fluid levels (alarm will sound if sheath tank is empty or waste is full)
 - o Sheath tank and plenum
 - o Waste tank
 - o Don't silence alarm
 - o Hit restart then hold prime button down until the light goes off
 - o Don't turn off FACFlow supply system
- Bleed out air in the bubble filter
- Instrument control
 - o Sample injection port (SIP)
 - Must use Falcon brand polystyrene tubes
 - Check for cracks
 - o Acquisition control panel
 - Standby and Run
 - Run means the sample is running regardless of what is going on in the software
 - Flow rate – don't use high button
 - Fine adjustment – keep between 250-350
 - Prime – avoid using unless you get a clog
 - Run di water after for 5 minutes if prime is used
 - Don't touch the mode button - OK if flashing

Instrument Setup

- Open FACSDiva software
 - o Lab login
 - o No password
 - o CST settings
 - Adjusts voltages and laser delays based on the daily morning QC
 - o Make a folder with your name
- Open experiment from template or create new
 - o Always rename experiments to include a date
- Instrument settings
 - o Detectors and voltages
 - FSC and SSC
 - Doublet discrimination
 - Channel names based on excitation laser and bandpass emission filter
 - Delete all fluorescent parameters and then add each one
 - Configuration on printout and ECFCC website
 - Voltrations
 - Stain index

Compensation controls

- Experiment menu -> create comp controls
 - o Comp controls specimen is created automatically
 - o Universal negative option
 - Provides a worksheet that displays all parameters at once
 - Run single stains first to make sure all are on scale and in correct detector
 - Bring beads or cells on scale in scatter plot
 - Apply gate to all compensation controls (if using all beads or all cells)
 - o May need to titrate Abs for the Symphony due to its high sensitivity
- Record compensation tubes
 - o Adjust scatter gate and positive gates – software should find the positive population but confirm gating and manually move gate to the peak fluorescence if necessary
 - o Keep universal negative or modify compensation to remove
- Calculate compensation
 - o Experiment menu -> calculate compensation
 - o Apply only

Experiment setup

- Adding specimens to experiment (syringe icon)
 - o Provides way of organizing tubes
 - o File naming

- Adding tubes
- Global vs Normal worksheets
- Add dot plots
 - Change parameters in plots to reflect fluorochromes being used
- Experiment layout
 - Create labels
 - Set acquisition stopping parameters
 - Acquire enough events to be statically significant
 - Can also set stopping parameters based on gates drawn (i.e. 100,000 singlets)

Acquisition

- Voltages
 - FSC and SSC voltages can be adjusted, especially between running beads and cells
 - Fluorescence voltages cannot be adjusted after compensation
- Acquire/Record
 - Restart button
 - Refresh rate of events on worksheet
 - Append vs overwrite and associated files
 - Gating, changing gate color
 - Show population hierarchy
 - Probe backflush between samples

Cleaning/shutdown

- Creating cleaning specimen from template
 - Run on high and record tubes
 - Do not overfill cleaning tubes with fluid
 - Okay if there are events in cleaning solution (doesn't mean it's dirty)
 - Instrument in standby and water on probe when finished cleaning
 - Shutdown if last user of the day

Exporting data

- Data automatically exported to PC
 - Desktop -> FCS folder
 - Transfer to OneDrive while cleaning
 - No flashdrives
 - Data stays on PC for one month
 - Experiments stay on PC for one month – create templates

Troubleshooting

- Fluid level isn't moving in tube
- Cytometer isn't connected
 - o Power cycle instrument and computer
 - o Windows password
 - o Turn UV laser back on

Resources and PPMS

- Symphony resources on ECFCC website
- Booking through PPMS
 - o 24 hour cancellation policy
 - o Rescheduling/deleting time
- Incident reporting
- Requesting assistance